


PDS No. 227261/ 227270 / 227280 227281 / 227283 227285	PRODUCT DATA SHEET				Page 1 of 1
Revision 04	Tube, 50 ml, PP, With Screw Cap				
	Item-No. 227261 / 227270 / 227280 / 227281 / 227283 / 227285				
Valid for Item-No.:	227261 (sterile)	227270 (sterile)	227280 (sterile)	227281 (sterile)	227283 (sterile)
	227285 (sterile)				

1.	Description / Specification	
1.1	Description	PP tube with screw cap, conical bottom, graduation and writing area 227261, -270, -280, -283: blue screw cap 227281, -285: white screw cap
1.2	Dimensions	Tube: 29,5 x 114,5 mm Tube with screw cap: 34,9 x 115,5 mm
1.3	Volume	Total volume: 50 ml
1.4	Material / Resin	Tube: PP (Polypropylene), free from heavy metal Cap: HDPE (High Density Polyethylene), free from heavy metal
1.5	Colour	Tube: 227261, -270, -281, -285: natural ; print: blue graduation, white writing area 227280, -283: brown ; print: graduation blue, writing area white Cap: 227261, -270, -280, -283: blue 227281, -285: white
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free from DNase/RNase and human DNA, non-pyrogenic, non-cytotoxic
2.2	Temperature range	Tube: -196°C to +121°C Cap: -50°C to +100°C
2.3	Autoclavability	Yes (without cap)
2.4	Centrifugation, max. RCF	3200 x g: Swinging-bucket rotor 9500 x g: Fixed-angle rotor
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging / Labelling	227261, -280, -285	227270, -281, -283
3.1	Pieces / Bag or Rack	20	25 in rack
3.2	Pieces / Box	500	300 pieces (12 racks)
3.3	Lot-No.	YYWWXXPP (year, week, index, working place or tool)	YYWWXXPP (year, week, index, working place or tool)
3.4	Other Information	QA-Certificate	QA-Certificate

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 15 August 2006	Date 15 August 2006	Date 15 August 2006	
Date 09.03.2005	Name T. Brando	Name Dr. C. Baumstark	Name Dr. C. Baumstark	