





96 Well ELISA Microplates U-Bottom / V-Bottom / F-Bottom Half Area

- Manufactured from crystal clear polystyrene

Free of detectable
DNase, RNase,
human DNA
non-pyrogenic



Cat.-No.	650 001	650 061	651 001	651 061
Well format	96 well	96 well	96 well	96 well
Well profile	U-bottom	U-bottom	V-bottom	V-bottom
Binding	MICROLON® 200 med. binding	MICROLON® 600 high binding	MICROLON® 200 med. binding	MICROLON® 600 high binding
Colour	clear	clear	clear	clear
Lid	-	-	-	-
Quantity per bag/case	5/40	5/40	5/40	5/40



Cat.-No.	655 001	655 061	655 080	655 081	675 001	675 061
Well format	96 well	96 well	96 well	96 well	96 well	96 well
Well profile	F-bottom / ST	F-bottom / ST	F-bottom / chimney well	F-bottom / chimney well	half area	half area
Binding	MICROLON® 200 med. binding	MICROLON® 600 high binding	MICROLON® 200 med. binding	MICROLON® 600 high binding	MICROLON® 200 med. binding	MICROLON® 600 high binding
Colour	clear	clear	clear	clear	clear	clear
Lid	-	-	-	-	-	-
Quantity per bag/case	5/40	5/40	5/40	5/40	10/40	10/40

1 Cell/
Tissue Culture

2 HTS-
Microplates

3 Immunology/
HLA

4 Microbiology/
Bacteriology

5 Tubes/Multi-
Purpose Beakers

6 Liquid
Handling

7 Molecular
Biology

8 Protein
Crystallisation

9 Separation

10 Biochips/
Microfluidics

11 Cryo-
Technics

12 Lids/Sealers/
CapMats

13 Reaction Tubes/
Analyser Cups

14 Accessories