



Containers, specimen, LDPE



Azlon

- Containers with captive hinged cap suitable for many applications
- Lightweight and unbreakable in normal use
- Can withstand temperatures up to 80°C

Product Code	Capacity (ml)	O.D. x height mm	Pack qty.
BGG302	1	8 x 38	20
BGG304	2.5	12 x 40	20
BGG306	5	15 x 49	20
BGG308	10	28 x 32	10
BGG311	20	32 x 36	10
BGG312	25	25 x 78	10



Containers, specimen, polystyrene



Azlon

- Clear polystyrene containers with maximum sample visibility
- Polyethylene airtight snap-on lids for quick easy use

Product Code	Capacity (ml)	O.D. (mm)	Height mm	Pack qty.
BGG440P	10	26	28	10
BGG442P	20	32	33	10
BGG444P	30	36	33	10
BGG446P	50	36	53	10
BGG450P	100	51	60	10
BGG452P	150	65	53	10
BGG454P	200	65	88	10
BGG456P	250	65	87	5
BGG458P	500	100	100	5