

BEAKER

Beakers, moulded in Polypropylene, have excellent clarity and very good chemical resistance. These autoclavable Beakers have prominently raised graduations on the outer surface. Pouring a solution becomes very convenient because of its tapered spout.

Material : Polypropylene



CODE NO.	PARTICULARS	LEAST COUNT	PACKING
11101	50 ml.	5	12 Pcs.
11102	100 ml.	5	12 Pcs.
11104	250 ml.	10	12 Pcs.
11105	500 ml.	10	12 Pcs.
11106	1000 ml.	25	6 Pcs.
11107	2000 ml.	250	6 Pcs.
11108	3000ml	300	6
11109	5000ml	500	4