

18cm Cell Scraper



70-1180

18cm Cell Scraper

FEATURES

- *18cm
- *Non-pyrogenic
- *Polyethylene
- *Pivoting blade makes precise cutting and scraping movements
- *Appropriate for a variety of cell culture vessels
- *The uniquely designed handle and blade reduce cell damage during the harvesting process
- *Pre-sterilized and individually wrapped in peel-open packaging

25cm Cell Scraper



70-1250

25cm Cell Scraper

FEATURES

- *25cm
- *Non-pyrogenic
- *Polyethylene
- *Pivoting blade makes precise cutting and scraping movements
- *Appropriate for a variety of cell culture vessels
- *The uniquely designed handle and blade reduce cell damage during the harvesting process
- *Pre-sterilized and individually wrapped in peel-open packaging

Cell Lifter



70-2180

Cell Lifter

FEATURES

*Used with dishes and cell culture clusters for quick removal of cell layers

*Blade Length: 2.0cm; Handle Length: 16cm

*Gamma radiation sterilized and individually wrapped

*Disposable lifters quickly remove cell layers

*The one-piece design and chiseled blade provide improved control while minimizing damage during the cell removal process

*Non-pyrogenic

*Polyethylene