

Tweezers, PMP / POM



- High quality, autoclavable tweezers
- Available with pointed or rounded ends
- Fine tooth grips for secure handling of samples

Product Code	Material	Colour	Ends	Pack qty.
TWY260	PMP	White	Sharp	1
TWY262	PMP	White	Sharp	1
TWY264	POM	Yellow	Round	1
TWY266	POM	Yellow	Round	1



FORCEPS

SAMCO Spring Forceps

Glass filled nylon, autoclavable to 125°C

Code	Size mm
L760/01	125

