


## ISOFLEX LABORATORY HOSES



Autoclavable PVC hoses, transparent non-toxic; non ageing and non oxidizing. Supplied in 20 meter rolls packed in PE bags. The hoses is marked every meter for ease of use.

### Technical data

Colour: clear

Odour: none

Taste: none

Toxicity: non-toxic

Ageing: none

Oxidation: none

### Chemical resistance

Strong acids: good

Weak acids: excellent

Strong alkalis: good

Weak alkalis: excellent

Organic solvents: fair

## Material:

PVC NON-TOXIC, AUTOCLAVABLE

Art.	I.D. Ø mm	O.D. Ø mm	Thickness mm
1760	0,8	2,4	0,8
1761	2,0	4,0	1,0
1762	4,0	6,0	1,0
1763	5,0	8,0	1,5
1764	6,0	9,0	1,5
1765	7,0	11,0	2,0
1766	8,0	12,0	2,0
1767	10,0	14,0	2,0
1768	12,0	17,0	2,5