


Tubing clamp

For fine control of liquids and gases

These thermoplastic polyester clamps are designed for use with flexible rubber and plastic tubing

- For use with flexible rubber and plastic tubing
- Provides fine flow control of liquids and gases
- Easy one-handed adjustment: just roll serrated wheel along sloped floor
- Manufactured from thermoplastic polyester



Item Code

For tubing size

KT4.5

O.D. up to 4.5mm

KT6

O.D. up to 6mm

KT10

O.D. up to 10 mm

KT14

O.D. up to 14mm