



SMART SPATULAS™

Product code :120003

Size-related design to fit perfectly with standard 15 and 1.5 mL tubes.

- Color: standard blue
- Ideal for dispensing or sample collection of powders or soft solids
- Shatter resistant in liquid nitrogen and resistant to dilute acids and bases and can be used with hot liquids
- Anti-static, sterile and eco-friendly versions available
- Autoclavable and complies with FDA standards for food contact
- Disposable

Upc Code	Code	Description	Use With Tubes	Other(s)
120003	120003	Standard Blue	15 and 1.5 mL	<ul style="list-style-type: none"> • L x W: 210 x 7 mm • Spoon (Flare): 25 x 14 mm • Scoop or Pick (l x h): 30 mm x 3 mm • Qty: 300/pk
120004	120004	Standard Opaque	15 and 1.5 mL	<ul style="list-style-type: none"> • L x W: 210 x 7 mm • Spoon (Flare): 25 x 14 mm • Scoop or Pick (l x h): 30 mm x 3 mm • Qty: 300/pk

Upc Code	Code	Description	Use With Tubes	Other(s)
120005	120005	Micro Antistatic	0.2 and 0.5 mL	<ul style="list-style-type: none"> • L x W: 140 x 3.5 mm • Spoon (Flare): n/a • Scoop or Pick (l x h): 11 mm L tapered @ 15° (V-Tip) • Qty: 300/pk •
120007	120007	Standard Opaque Sterile	15 and 1.5 mL	<ul style="list-style-type: none"> • L x W: 210 x 7 mm • Spoon (Flare): 25 x 14 mm • Scoop or Pick (l x h): 30 mm x 3 mm • Qty: 100/pk
120006	120006	Macro Green	50 mL	<ul style="list-style-type: none"> • L x W: 310 x 10 mm • Spoon (Flare): 42 x 23.5 mm • Scoop or Pick (l x h): 30 mm L tapered @ 35° • Qty: 150/pk
120008	120008	Eco	15 and 1.5 mL	<ul style="list-style-type: none"> • L x W: 140 x 7 mm • Spoon (Flare): 25 x 14 mm • Scoop or Pick (l x h): 30 mm x 3 mm • Qty: 300/pk