



Rack, microcentrifuge tube, polypropylene *A E Azlon*

- Economical and sturdy test tube racks supplied in self assembly flat packs
- Accepts 128 x 1.5ml microcentrifuge tubes
- Autoclavable at 121°C
- Sinks in water baths
- Alpha numeric place coding

Product Code	Colour	l x w x h mm	No. places	Pack qty.
RPB001P	Blue	265 x 126 x 38	128	5
RPR001P	Red	265 x 126 x 38	128	5
RPW001P	White	265 x 126 x 38	128	5
RPY001P	Yellow	265 x 126 x 38	128	5



Racks, microcentrifuge tube, polypropylene *A E Azlon*

- Microcentrifuge tube rack with numerical marking
- Autoclavable at 121°C
- Holds a maximum of 20 tubes

Product Code	Base dim. l x w x h mm	Hole diam. mm	Pack qty.
RWK058	210 x 70 x 35	10	1

