


# GENERAL APPARATUS

## TUBING CLIPS



### Ratchet Clip

Closes tubing to external diameters shown

Code	Size mm	Material
L140/02	10	Polypropylene
L140/03	17	Polyester

### Hoffman Clip

Bunsen screw clip, nickel plated brass fitted with protection plate and hinged bottom bar. Dimensions shown are width of gate

Code	Size mm
L141/01	13
L141/02	20
L141/03	25
L141/04	32
L141/05	40

### Hoffman Clip

Heavy duty fitted with nickel plated brass tommy screw and hinged bottom bar. Dimension shown is width of gate

Code	Size mm
L142/02	40

### SAMCO Screw Tubing Clip

Closes rubber and flexible PVC tubing of up to 14mm external diameter. The hinged side allows it to be clamped around the tubing, already connected to the apparatus. The screw adjustment ensures fine flow control. Moulded in acetal copolymer which is easy to clean, can be sterilised and is resistant to most acids and alkalis

#### Code

L145/01

### SAMCO Spring Tubing Clip

Operates with a parallel action to provide optimum closure of rubber and flexible PVC tubing. Made from chemical and stain resistant acetal copolymer and stainless steel.

Size 8: accepts tubing up to 8mm outside diameter  
Size 12: accepts tubing up to 12mm outside diameter

#### Code

Size  
mm

L150/01

8

L150/02

12