

SciLabware Ltd., Unit 4, Riverside 2, Campbell Road, Stoke on Trent, ST4 4RJ, UK Tel; +44 (0)1782 444 406
Fax; +44 (0)1782 940 436 info@scilabware.com www.scilabware.com

2nd March 2015

Azlon® new polypropylene wash bottles

New products for 2015

SciLabware is pleased to introduce a new range of Azlon polypropylene, non-vented, wash bottles. The ability of polypropylene to withstand higher temperatures and resist some solvents better than LDPE, make Azlon polypropylene wash bottles the obvious choice for use in the following applications:

Autoclavable wash bottle

The Azlon autoclavable wash bottle is manufactured entirely from polypropylene and is suitable for use with autoclaves operating at 121°C.

Dichloromethane wash bottle

Our new multi-lingual, non-vented, Dichloromethane wash bottle is printed with the chemical name, CAS number, NFPA diamond, Globally Harmonised System (GHS) symbols and safety information in the following four languages, English, French, German and Spanish.

Methyl-ethyl-ketone (MEK) wash bottle

Our new multi-lingual, non-vented, Methyl-ethyl-ketone wash bottle is printed with the chemical name, CAS number, NFPA diamond, GHS symbols and safety information in the following four languages, English, French, German and Spanish.

Each of these products is available in our popular 500ml, easy fill, wide mouth design.

One of our sales representatives will contact you to assist you with these and other products from our range.

If you require any further information for these or any other SciLabware products then please visit www.scilabware.com

Kind regards



Michael Fowler
Product Manager – Azlon Plastics
michael.fowler@scilabware.com



Azlon® new polypropylene wash bottles



Bottles, wash, autoclavable, non-venting, wide neck, PP



- Autoclavable at 121°C
- PP bottle fitted with one-piece PP closure
- Closures feature precision jet-tip aperture for fine control

Product Code	Capacity (ml)	Print label	Closure colour	Pack qty.
WGW700PP	500	Autoclavable	White	5



Caution. Do not autoclave if containing solvents or hazardous/reactive solutions or substances.



Caution. When autoclaving this item always ensure the cap is completely loosened or removed to prevent accidental implosion.



Bottles, wash, multi-lingual, non-venting, wide neck, PP



- Multi-lingual wide neck wash bottles
- Indelibly printed in four languages (English, French, German and Spanish)
- Solvent labels including the new GHS symbols, NFPA diamond, risk phrases and safety information
- Fitted with colour coded one-piece PP closures that maintain their shape in use
- Closures feature precision jet-tip aperture for fine control

Product Code	Capacity (ml)	Print label	Closure colour	Pack qty.
WGW701PP	500	Dichloromethane	White	5
WGW703PP	500	Methyl ethyl ketone	Orange	5



Caution. Do not autoclave WGW701PP and WGW703PP.

