

Wash bottle, LDPE, Bleach



- Multi-lingual wide neck wash bottle
- Opaque white LDPE for light sensitive sodium hypochlorite / bleach
- Indelibly printed in three languages - English, French and Spanish
- Features the new Globally Harmonised Symbols, NFPA diamond, risk phrases and safety information
- Fitted with colour coded one-piece yellow PP closure
- Closure features precision jet-tip aperture for fine control

Item Code

Capacity ml

Printing label

WG401PML

500

Sodium Hypochlorite