



Ramboldi® High Temperature Autoclave Bags

Autoclave bags for **hazardous waste**

Ramboldi autoclave bags are available in three different sizes and are printed with blue hazard warnings labels for contaminated waste disposal

- Ideal for the disposal of all biological waste in autoclaves or incinerators
- Manufactured from strong polypropylene
- Suitable for use up to 135°C
- Uniform wall thickness and seamless bottom, designed to ensure high pressure bearing
- Printed with blue bio-hazard labelling specifically designed for contaminated waste disposal
- Write on panel to enable the user to mark each bag as required



Product Specification

Item Code	Bag Material	Description	Inner Pack Qty	Label	Length	Width mm	Width x Length (mm)	Pack quantity
370100	HDPE	High Temperature Autoclave Bag, Printed, 310 x 660mm	50	Printed	660mm	310	310 x 660	200
370110	HDPE	High Temperature Autoclave Bag, Printed, 415 x 600mm	50	Printed	600mm	415	415 x 600	200
370120	HDPE	High Temperature Autoclave Bag, Printed, 610 x 810mm	50	Printed	810mm	610	610 x 810	200
350190		Autoclave tape			50m	19		18
350195		Autoclave tape			50m	25		36