



### Multi-Channel Reagent Reservoirs

#### Designed and packaged for superior handling!

Our redesigned reservoirs are improved in every way \_ with thicker walls, crush-resistant "turrets" at the corners, improved pour spouts, 25 mL graduation marks, and tabs to prevent stacked reservoirs from sticking together. The bottom of each reservoir tapers down into a narrow trough, to help you get the last drop of liquid.

Sterile reservoirs are sold in packs of 50. Individually sealed in our new, Sani-Tear, high density polyethylene bag and gamma-irradiated in an ISO-certified US facility.

Non-sterile reservoirs are sold in packs of 100.

<b>20521A</b>	55 mL	Clear PVC	Non-Sterile	Case of 16		
				Pckg. of 100		
<b>20521B</b>	55 mL	White Polystyrene	Non-Sterile	Case of 16		
				Pckg. of 100		
<b>20521C</b>	55 mL	White Polystyrene	Sterile	Case of 6		
				Pckg. of 50		
<b>20721B</b>	100 mL	White Polystyrene	Non-Sterile	Pckg. of 100		
<b>20721C</b>			Sterile	Case of 6		
				Pckg. of 50		