

Low form Griffin style beakers, PMP



- Complies with ISO 7056
- Traditional beaker shape
- Ideal for general laboratory use
- Suitable for food contact applications in accordance with U.S. FDA regulations
- Good chemical resistance

- Can withstand temperatures up to 150°C and therefore autoclavable
- Graduations are clearly printed in indelible solvent resistant ink

Item Code	Capacity ml	Diameter mm	Graduation ml	Height mm	Pack quantity
BDD312P	25	35	5	50	12
BDD318P	50	42	10	60	12
BDD324P	100	51	10	69	12
BDD330P	250	71	25	94	6
BDD334P	400	80	50	111	6
BDD336P	500	87	50	117	6
BDD338P	600	95	50	125	6
BDD340P	1000	107	100	148	6
BDD342P	2000	131	200	187	6