



**i** Exposure to temperatures above 80°C may affect accuracy.

## Cylinders, moulded graduations, class B, polypropylene



Azlon

- Conform to ISO 6706
- Ultra clear polypropylene with good chemical resistance
- Large pouring spout for easy dispensing
- Stable hexagonal bases with non slip anti-suction feature
- Permanent moulded graduations

Product Code	Capacity (ml)	Grads. ml	Tol. $\pm$ ml	Height mm	Pack qty.
CP0010P	10	0.2	0.20	140	10
CP0025P	25	0.5	0.50	170	10
CP0050P	50	1.0	1.00	200	10
CP0100P	100	1.0	1.00	250	5
CP0250P	250	2.0	2.00	215	1
CP0500P	500	5.0	5.00	360	1
CP1000P	1000	10.0	10.00	440	1
CP2000P	2000	20.0	20.00	535	1

## Cylinders, moulded graduations, set, class B, polypropylene

Product Code	Description	Pack qty.
NB1203	Set of 6 cylinders comprising (1x 10, 25, 50, 100, 250 & 500ml)	1



## Cylinders, printed graduations, class A, PMP (TPX)



Azlon

- Calibrated to ISO 6706 and DIN 12681, class A tolerances
- Steam sterilisation at 121°C does not affect accuracy
- Supplied with conformity certificate
- Enhanced pin sharp blue printed graduations
- Stable hexagonal bases with non-slip anti-suction feature

Product Code	Capacity (ml)	Grads. ml	Tol. +/- ml	Height mm	Pack qty.
CTA0010P	10	0.2	0.1	140	2
CTA0025P	25	0.5	0.2	170	2
CTA0050P	50	1.0	0.5	200	2
CTA0100P	100	1.0	0.5	250	2
CTA0250P	250	2.0	1	315	2
CTA0500P	500	5.0	2.5	360	2
CTA1000P	1000	10.0	5	440	1
CTA2000P	2000	20.0	1	535	1