



Cylinders, printed graduations, class B, PMP (TPX)



- Calibrated to ISO 6706 and DIN 12681, class B tolerances
- Steam sterilisation at 121°C does not affect accuracy
- Enhanced pin sharp blue printed graduations
- Stable hexagonal bases with non-slip anti-suction feature

Product Code	Capacity (ml)	Grads. ml	Tol. +/- ml	Height mm	Pack qty.
CTB0010P	10	0.2	0.2	140	10
CTB0025P	25	0.5	0.5	170	5
CTB0050P	50	1.0	1	200	5
CTB0100P	100	1.0	1	250	5
CTB0250P	250	2.0	2	315	1
CTB0500P	500	5.0	5	360	1
CTB1000P	1000	10.0	10	440	1
CTB2000P	2000	20.0	20	535	1



Cylinders, moulded graduations, class B, PMP (TPX)



- Conform to ISO 6706
- Transparent robust PMP
- Good chemical resistance
- Large pouring spout for easy dispensing
- Autoclavable to 121°C

Product Code	Capacity (ml)	Grads. ml	Tol. +/- ml	Height mm	Pack qty.
CT0010P	10	0.2	0.2	140	10
CT0025P	25	0.5	0.5	170	5
CT0050P	50	1.0	1	200	5
CT0100P	100	1.0	1	250	5
CT0250P	250	2.0	2	315	1
CT0500P	500	5.0	5	360	1
CT1000P	1000	10.0	10	440	1
CT2000P	2000	20.0	20	535	1

