



Racks, draining, coated




Azlon

- High quality free standing or wall mounted draining racks
- Manufactured from polyamide epoxy coated steel
- Offers a better heat and chemical resistance than other coatings
- There are eight pegs 110mm long and twenty four pegs 80mm in length. Peg diameter is 4.5mm
- Overall dimensions of the rack are 410 x 300mm

Product Code	Description	Pack qty.
RWK070FS	Rack with integral free standing support	1
RWK070	Wall mountable	1

 *RWK070 shown.*



Wall mountable draining rack, polystyrene



Azlon

- Lightweight, easy to assemble polystyrene draining racks designed to be wall mounted
- Pegs simply slot into place on the back board and are spaced to allow various sizes of labware to fit on to the rack
- Water is collected in the trough at the base of the rack and drains away via a rubber tube (not supplied) which can be directed into a nearby sink or container
- Supplied with wall mounting accessories
- Measures approximately 650 mm high x 450 mm wide
- A set of narrow pegs for test tubes, volumetric flasks etc. is also available

Product Code	Description	Peg l x O.D. mm	Pack qty.
WRP001P	Rack, including 72 pegs, trough & drain	120 x 16	1
WRP002P	Set of 11 narrow pegs	100 x 6	1

