



Flasks, conical, vacuum



- With heavy walls and bases
- Suitable for use under vacuum

Product Code	Capacity (ml)	Socket size	Approx height (mm)	Approx diam. mm	Pack qty.
FEH250/3	250	24/29	139	90	1



Always take safety precautions when using flasks under vacuum.




Flask, conical, polypropylene



- Polypropylene erlenmeyer flask with excellent chemical resistance
- Wide neck with PP screw cap
- Flask and cap can be autoclaved at 121°C

Product Code	Capacity (ml)	Grads. ml	Thread size	Neck size	Pack qty.
FDE542	50	10	GL40	34/35	1
FDE544	100	20	GL40	34/35	1
FDE546	250	50	GL52	45/40	1
FDE548	500	100	GL52	45/40	1
FDE550	1000	200	GL52	45/40	1



Flask, conical, PMP



- Manufactured from transparent PMP
- Wide neck with PP screw cap
- Flask and cap can be autoclaved at 121°C

Product Code	Capacity (ml)	Grads. ml	Thread size	Neck size	Pack qty.
FDE442	50	10	GL40	34/35	1
FDE444	100	20	GL40	34/35	1
FDE446	250	50	GL52	45/40	1
FDE448	500	100	GL52	45/40	1
FDE450	1000	200	GL52	45/40	1