

Polypropylene Vacuum Flask



Related applications: [General Plasticware](#)

Replace glass flasks with the Polypropylene Vacuum Flask, molded in one piece for full vacuum service with angled tubulation to prevent tipping. This shatter-proof, autoclavable flask can withstand vacuum to 30" Hg

Catalog number	Product Size	Capacity (Metric)	Accepts
DS4101-0500	Each	500mL	#7 stopper
DS4101-1000	Each	1000mL	#8 stopper
DS4101-2000	Each	1700mL	#10 stopper