



Ice bucket with lid, polyurethane



Azlon

- Black, green and blue skinned polyurethane foam ice buckets
- Suitable for use with water-ice, dry-ice and salt solutions
- Ideal for separating work areas by colour
- Supplied with lids
- Integral handles allow easy lifting and carrying
- Octagonal shape for easy stacking and stability
- Capacity 4.5 litres
- Dimensions 185 x 330mm (h x w) and weight 1200g

| Product Code | Colour | Pack qty. |
|--------------|-----------|-----------|
| IBB001P | Black | 1 |
| IBB002P | Pale Blue | 1 |
| IBB003P | Green | 1 |



Not suitable for use with solvents or liquid nitrogen.



Mini Ice Bucket

Azlon

- Skinned polyurethane foam mini ice buckets
- Suitable for use with water-ice, dry-ice and salt solutions
- Integral handles allow easy lifting and carrying
- Tapered for easy stacking and stability
- Capacity 1000ml
- Dimensions 95mm x 200mm (h x w) and weight 140g

| Product Code | Colour | Pack qty. |
|--------------|--------|-----------|
| IBB004P | Green | 1 |



Not suitable for use with solvents or liquid nitrogen.