

WIDE MOUTH STORAGE BOTTLES

Material: High density Polyethylene

Color: Neutral- For food use

They are equipped with inner-cap and cap with holes for warranty seal.

Very strong and durable for extended use over time. They are equipped with inner-cap for a tight seal. Supplied with cap with holes for warranty seal. Ideal for containing liquids, powders and granulates.

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
PL 416	15	240	460	115
PL 417	20	265	500	130
PL 1350	25	275	570	130

