

# Jugs

## Jugs, printed graduations, polypropylene *A E Azlon*

- Tapered jugs in high clarity polypropylene
- Suitable for food contact applications in accordance with EU and U.S. FDA regulations
- Crisp blue printed graduations
- Autoclavable at 121°C
- Ergonomically designed handle

Product Code	Capacity (ml)	Rim diam. (mm)	Height (mm)	Pack qty.
JPB0050P	50	45.0	68	5
JPB0100P	100	61.0	72	5
JPB0250P	250	70.0	120	1
JPB0500P	500	90.0	140	1
JPB1000P	1000	104.0	180	1
JPB2000P	2000	150.0	214	1
JPB3000P	3000	152.0	240	1
JPB5000P	5000	200.0	270	1



## Jugs, moulded graduations, polypropylene *A E Azlon*

- Tapered jugs in high clarity polypropylene
- Suitable for food contact applications in accordance with EU and U.S. FDA regulations
- Permanent moulded graduations
- Autoclavable at 121°C
- Ergonomically designed handle

Product Code	Capacity (ml)	Sub div. (ml)	Rim diam. (mm)	Height (mm)	Pack qty.
JPM0050P	50	2	45.0	68	5
JPM0100P	100	2	61.0	72	5
JPM0250P	250	5	70.0	120	1
JPM0500P	500	10	90.0	140	1
JPM1000P	1000	10	104.0	180	1
JPM2000P	2000	50	150.0	214	1
JPM3000P	3000	50	152.0	240	1
JPM5000P	5000	100	200.0	270	1

