

AZLON® SCREW CAP CONTAINERS, POLYPROPYLENE

Perfect for the **storage** of powders and samples

Azlon® screw cap containers are composed of polypropylene which provides good chemical resistance and are available in sizes from 30ml to 1.2L.

- Wide neck containers with liner-less screw caps
- Good chemical resistance
- Ideal for storage of samples, powders, etc.



Product Specification

Item Code	Diameter mm	Capacity ml	Height mm	Pack quantity
JWP0030P	43	30	43	72
JWP0060P	53	60	47	48
JWP0125P	70	125	60	36
JWP0250P	89	250	66	36
JWP0500P	89	500	97	24
JWP0950P	120	950	100	24
JWP1200P	120	1200	132	24