



Lead Rings

Use rings to stabilize glassware and plasticware on shakers or immersed in water baths. Vinyl-coated lead rings will not mark or scratch bench tops or labware.

Cat No.	Type	Flask Sizes	ID mm
HD8882A	O-Shape	125 to 500 mL	48
HD8882B	O-Shape	250 to 1000 mL	51
HD8882C	O-Shape	500 to 2000 mL	57
HD8882D	O-Shape	1000 to 4000 mL	70
HD8881A	C-Shape	125 to 500 mL	42
HD8881B	C-Shape	250 to 1000 mL	51
HD8881C	C-Shape	500 to 2000 mL	54
HD8881D	C-Shape	1000 to 4000 mL	66