



## 100-WELL MICROTUBE STORAGE BOX, ASSORTD, 5 / PK

Product Code:120044

UPC Code: 120044

100-Well Microtube Storage Boxes

Secure storage for microtube and cryogenic vials

- For larger storage requirements of 1.5/2.0 mL microtubes or cryogenic vials.
- Black, light-tight version available for light sensitive samples
- Easy-open hinged lid with snap clasp
- 100 wells
- Well dia. 1.23 cm
- Dimensions: 5.6 x 6 x 2.2 in (14.1 x 15.1 x 5.7 cm)
- 5/pk and 8/cs
- Color: Assorted (blue, green, purple, yellow, orange)
- Autoclavable



## 50 WELL MICROTUBE STORAGE BOX

Product Code:120033

UPC Code:HD120033

Ideal for sample or case study storage.

- Made from polypropylene
- Designed for quick visual identification with uniform wells and location identification
- Tube storage capacity is indicated on the outside of the box
- Translucent vibrant colors with a frosted labeling area
- Debossed with a distinctive wave design
- Stackable, saves valuable space in the lab
- Alphanumeric Identification for quick location and retrieval of tubes
- Store down to  $-80^{\circ}\text{C}$ . - Must use tubes capable of being stored at the required temperature of the user
- Great for smaller storage requirements of 1.5/2.0 mL microtubes or cryogenic vials
- Easy-open hinged lid with snap clasp
- Tubes: 1.5/2.0 mL
- Wells: 50
- Well Dia: 12.3 mm
- Dimensions: 5.6 x 3.6 x 2.2 in (14.2 x 9.2 x 5.6 cm)
- Packaging: 5ea/pk and 20pk/cs
- Autoclavable



## 81-WELL MICROTUBE STORAGE BOX,

Product Code:120040

---

120040      120040      Assorted (Blue, Green, Purple, Yellow, Orange)

### 81-Well Microtube Storage Boxes

Compact footprint for microtube and cryogenic vial storage.

- Made from polypropylene
- Store 1.5/2.0 mL microtubes or cryogenic vials in polypropylene storage boxes
- Easy open, friction fit lid
- Designed for quick visual identification with uniform wells and location identification
- Tube capacity is indicated on the outside of the box
- Translucent vibrant colors with a frosted labeling area
- Alphanumeric Identification for quick location and retrieval of tubes
- Stackable, saves valuable space in the lab
- Debossed with a distinctive wave design
- Store down to -80°C.- Must use tubes capable of being stored at the required temperature of the user
- Black light tight version available for light sensitive samples
- Tubes: 1.5/2.0 mL
- Wells: 81
- Well Dia: 12.3 mm
- Pkg: 5/pk and 4/cs
- Dimensions: 5.1 x 5.1 x 1.9 in (13 x 13 x 4.7 cm)
- Autoclavable