



Stand, petri dish, coated



- Holds up to six 90mm petri dishes
- Manufactured from polyamide epoxy coated steel
- Offers a better heat and chemical resistance than other coatings
- Ideal for small batch incubation

Product Code	Pack qty.
SWP062	1



Stand, retort, weighted, polyethylene



- Weighs 2.5kg ensuring stability
- Corrosion resistant
- Polyethylene coated metal rod screws into the base
- Available with the hole in the centre or at the short edge

Product Code	w x l x h mm	Rod position	Rod size mm	Pack qty.
SWP012	171 x 248 x 52	Centre	760 x 12.5	1
SWP014	171 x 248 x 52	Edge	760 x 12.5	1



i A range of burette clamps are featured on page 88.