



Pipette rack, pipette stand



Pipette rack, HDPE



Azlon

- White moulded polyethylene rack
- Will support 12 pipettes over 200 mm long
- Very stable
- Supplied as a kit of two end pieces and three connecting rods for easy push fit assembly

Product Code	Pack qty.
SWP038	1



Pipette stand, polypropylene



Azlon

- Blue chemically resistant polypropylene stand.
- Holds up to 94 pipettes vertically
- Very stable with 230mm diameter base.
- Height of 470mm
- Has a rotatable rack suspended by a support rod

Product Code	Description	Pack qty.
SWP033	Pipette Stand	1



Pipette support, polystyrene



Azlon

- Moulded in white polystyrene
- Simple and efficient pipette rest
- An off-centre support means it is also suitable for short pipettes, burettes, syringes etc.

Product Code	l x w x h mm	Pack qty.
SWP036	280 x 215 x 38	1