



Cube Rack

Interlocking on all sides!

Our user friendly cube rack is a must for any laboratory. Slide these polypropylene racks together in multiple orientations—work with different sized tubes at the same time. Accommodates (4) 50mL tubes, (10) 15mL conical tubes, (12) 12 x 75 or 12 x 100 tubes, or (16) 1.5 to 2.0mL microcentrifuge tubes. Autoclavable. Available in natural, blue, pink, yellow, or assorted colors (one each of blue, green, yellow, orange, and pink). Rack measures 4 1/4 x 4 1/4 x 4 1/2". Packs of 5

HD29050A	Cube Rack, Assorted Colors
-----------------	----------------------------



4-Way Microtube Rack

Link together multiple racks

Create as large a rack as you need by linking together multiple racks. Heavy-duty polypropylene racks feature a unique system of tabs and slots, which facilitate easy connection and sturdy fit. Each rack can hold four 50 mL conical tubes, twelve 15 mL conical tubes, thirty-two 1.5 mL microtubes, or thirty-two 0.5 mL microtubes. Autoclavable. Assorted colors: blue, green, pink, yellow, and orange. Rack measures 174 x 95 x 52 mm. Pack of 5.

Cat No.	Description	Units
HD29022A	4-Way Microtube Rack, Natural	Case of 20
		Pckg. of 5
HD29022B	4-Way Micortube Rack, Blue	Case of 20
		Pckg. of 5
HD29022C	4-Way Micortube Rack, Green	Case of 20
		Pckg. of 5
HD29022D	4-Way Microtube Rack, Pink	Case of 20