



80-Well Microtube Rack

Popular fraction-collector-plate style racks

Easily fill and transfer 80 tubes in compact arrangement (5 x 16 array). Polypropylene rack holds 1.5 to 2.0 mL microtubes. Molded-in coordinates aid organization. Autoclavable. Assorted colors: blue, green, pink, yellow, and orange. Measures 225 x 67 x 28 mm. Pack of 5

Cat No.	Description	Units
HD29025A	80 Well Microtube Rack, Natural	Case of 80
		Pckg. of 5
HD29025B	80-Well Microtube Rack, Blue	Case of 80
		Pckg. of 5
HD29025C	80-Well Microtube Rack, Green	Case of 80
		Pckg. of 5
HD29025D	80-Well Microtube Rack, Pink	Case of 80
		Pckg. of 5
HD29025E	80-Well Microtube Rack, Yellow	Case of 80
		Pckg. of 5
HD29025F	80-Well Microtube Rack, Orange	Case of 80
		Pckg. of 5
HD29025G	80-Well Microtube Rack, Assorted Colors	Case of 80
		Pckg. of 5