



## 96-Well Reversible Rack

### Easy-to-grasp tabs

Polypropylene rack has 96 wells for 0.5 mL microcentrifuge tubes on one side and 96 wells for 1.5 to 2.0 mL microcentrifuge tubes on the other side. Wells are arranged in standard 8 x 12 array. Handy lid snaps into place via tabs on each end of rack. 0.5 mL wells are separated by molded-in grid lines, making them easier to see, and raised rings are around the 1.5 to 2.0 mL wells make tubes easier to remove. Autoclavable. Assorted colors: blue, green, pink, yellow, and orange. Rack measures 246 x 121 x 50 mm. Pack of 5.

Cat No.	Description	Units
<b>HD2345A</b>	96-Well Reversible Racks, Assorted Colors	Case of 25
		Pckg. of 5
<b>HD2345B</b>	96-Well Reversible Racks, Natural	Case of 25
		Pckg. of 5
<b>HD2345C</b>	96-Well Reversible Racks, Blue	Case of 25
		Pckg. of 5
<b>HD2345D</b>	96-Well Reversible Racks, Green	Case of 5
		Pckg. of 5
<b>HD2345E</b>	96-Well Reversible Racks, Pink	Case of 5
		Pckg. of 5
<b>HD2345F</b>	96-Well Reversible Racks, Yellow	Case of 5
		Pckg. of 5