



Dual Solution™ Reservoir

Two basins in one. Innovative shape guides and protect pipettor tips

Polypropylene solution basin/reagent reservoir has a standard 50 mL basin on one side and 12 individual 5 mL basins on the other side. Wells are angled both downward and outward, guiding multichannel pipettor tips toward a centered position, which helps to avoid contamination and waste. Imprinted indices identify samples; loose-fitting lids prevent evaporation and contamination. Non-sterile. Pack of 25 reservoirs and 25 lids.

HD20821A	Dual Solution Reservoir	Pckg. of 25
-----------------	-------------------------	-------------