



## Scoops, measuring, polypropylene



- Made from food safe polypropylene
- Top is flat for levelling off powders, base is flat enabling scoops to stand horizontally
- The capacity is clearly moulded on each size

Product Code	Capacity (ml)	Overall length mm	Pack qty.
SWL002	2	60	5
SWL006	5	82	5
SWL008	10	102	5
SWL012	25	137	5
SWL018	50	160	5
SWL024	100	198	5
SWL030	250	259	5
SWL036	500	315	5
SWL140	1000	385	1



## Scoops, HDPE



- Sturdy scoops in natural rigid polyethylene
- Lightweight with good chemical resistance
- Food safe

Product Code	Max. capacity ml	Overall length mm	Pack qty.
SWL028	150	250	5
SWL033	1000	350	5
SWL040	1500	405	1