



SHARPSGUARD® 0.5 web

Our SHARPSGUARD® 0.5 with web opening, is ideal for community healthcare workers and easily fitted into a bag or glove compartment. This container is supplied with lid assembled.

- **Usable Capacity:** 0.4 Litres
- **Container weight (g):** 88.0
- **Container Size:** H206mm x L79mm x W45mm
- **Aperture Dimensions:** W26mm x L30mm
- **Carton contains:** 60
- **Carton Dimensions:** H288mm x L687mm x W368mm
- **Weight of carton and containers (kg):** 6.03



SHARPSGUARD® 2.5

Our SHARPSGUARD® 2.5 is an ideal general-purpose sharps container, it provides maximum capacity in the minimum working space. The SHARPSGUARD® orange 2.5 also fits into our SHARPSGUARD® tray range.

- **Usable Capacity:** 2.21 Litres
- **Container weight (g):** 232.0
- **Container Size:** H166mm x L176mm dia
- **Aperture Dimensions:** W77mm x L46mm
- **Carton contains:** 48
- **Carton Dimensions:** H346mm x L663mm x W343mm
- **Weight of carton and containers (kg):** 11.68



SHARPSGUARD® 5

Our SHARPSGUARD® 5 is an ideal general-purpose sharps container, it provides maximum capacity in the minimum working space.

- **Usable Capacity:** 3.98 Litres
- **Container weight (g):** 338.0
- **Container Size:** H285mm x L176mm dia
- **Aperture Dimensions:** W77mm x L46mm
- **Carton contains:** 48
- **Carton Dimensions:** H346mm x L1080mm x W358mm
- **Weight of carton and containers (kg):** 15.71



SHARPSGUARD® 11.5

Our SHARPSGUARD® 11.5 is a multi-purpose container featuring a robust pail handle for easy transportation. A wider aperture allows easy insertion of larger items.

- **Usable Capacity:** 9.27 Litres
- **Container weight (g):** 622.0
- **Container Size:** H249mm x L275mm dia
- **Aperture Dimensions:** W148mm x L48mm
- **Carton contains:** 20
- **Carton Dimensions:** H566mm x L943mm x W283mm
- **Weight of carton and containers (kg):** 13.23



SHARPSGUARD® 22

Our SHARPSGUARD® 22 is ideal for theatres, laboratories and sluice rooms with high volumes of hazardous waste. The wide aperture allows economical and convenient disposal of large items.

- **Usable Capacity:** 17.86 Litres
- **Container weight (g):** 943.0
- **Container Size:** H343mm x L327mm dia
- **Aperture Dimensions:** W148mm x L58mm
- **Carton contains:** 10
- **Carton Dimensions:** H326mm x L1148mm x W323mm
- **Weight of carton and containers (kg):** 10.21