



Made of polypropylene

This handy autoclavable rack can hold up to **50** cryogenic vials. Now with one hand, you can easily unscrew a Simport Cryovial® closure.

Thanks to an innovative universal locking system, the vials will securely lock in each well and will not turn.

Each position is identified with an alphanumeric index. Strong handles make it easy and safe to carry.

It is supported by five anti-skid rubber feet.

The rack is compact and stackable.

- **RED Color**
- Also available in three colors.
- Size: 10 cm x 20 cm x 25 mm H (4 x 8 x 1 in. H)