



PRODUCT DESCRIPTION

STOP COCKS

These polypropylene stop cocks are fitted with a PTFE plug & have serrated tubulations the accept from 6mm to 10mm tubing. These stop cocks are absolutely leak proof and are available in two different options i.c. 2 way & 3 way.

Material : Polypropylene / PTFE

| CODE NO. | PARTICULARS | PACKING (Pcs) |
|----------|------------------|---------------|
| 55301 | Stop Cocks 2 Way | 6 |
| 55302 | Stop Cocks 3 Way | 6 |