

Silicone Rubber (VMQ), meets the following standards:

BFR XV, FDA CFR 21§ 177.2600, Europ.Pharmacopoeia

Test for biocompatibility according to 10993-1

transparent, Durometer 60 Shore A

temperature range -60°C to +200°C, short periods up to +260°C

Inner-Ø mm	Outer-Ø mm	Wall Thickness mm	Roll Length m	
0,5	1,3	0,4	25	
0,5	2,5	1	25	
0,8	1,6	0,4	25	
1	1,8	0,4	25	
1	3	1	25	
1	4	1,5	25	
1	6	2,5	25	
1,5	2,3	0,4	25	
1,5	3	0,75	25	
1,5	3,5	1	25	
1,5	4	1,25	25	
2	2,6	0,3	25	
2	2,8	0,4	25	
2	4	1	25	
2	5	1,5	25	
2	6	2	25	
2,5	3,3	0,4	25	
2,5	4,5	1	25	
3	3,8	0,4	25	
3	5	1	25	
3	6	1,5	25	
Inner-Ø mm	Outer-Ø mm	Wall Thickness mm	Roll Length m	

3	7	2	25	
4	5	0,5	25	
4	6	1	25	
4	7	1,5	25	
4	8	2	25	
4	9	2,5	25	
5	6	0,5	25	
5	7	1	25	
5	8	1,5	25	
5	9	2	25	
5	10	2,5	25	
6	8	1	25	
6	9	1,5	25	
6	10	2	25	
7	9	1	25	
7	10	1,5	25	
7	11	2	25	
8	10	1	25	
8	11	1,5	25	
Inner-Ø mm	Outer-Ø mm	Wall Thickness mm	Roll Length m	
8	12	2	25	
8	14	3	25	
8	15	3,5	25	
9	12	1,5	25	
9	13	2	25	
10	13	1,5	25	
10	14	2	25	
10	15	2,5	25	

10	16	3	25	
10	18	4	25	
11	15	2	25	
12	16	2	25	
12	18	3	25	
13	19	3	25	
14	20	3	25	
15	21	3	20	
16	21	2,5	20	
18	24	3	20	
20	26	3	20	
25	33	4	15	
Vacuum Quality				
Inner-Ø mm	Outer-Ø mm	Wall Thickness mm	Roll Length m	
3	9	3	25	
4	12	4	25	
5	15	5	15	
6	12	3	25	
7	13	3	25	
7	15	4	25	
7	17	5	15	
8	16	4	25	
8	18	5	15	
10	20	5	15	
12	22	5	15	