



Bottles, round, wide neck, PTFE



- Opaque PTFE bottles
- Rigid and non-porous with a smooth surface
- Resist all chemicals except molten alkali metals and fluorine
- Slope shoulder for easy pouring of powders or liquids
- Can withstand temperatures upto 260°C

Product Code	Capacity (ml)	Diam. x height mm	Neck I.D. mm	Pack qty.
BWN118	50	43 x 76	25	1
BWN124	100	52 x 88	35	1
BWN130	250	67 x 115	42	1



Do not use when sealed for chemical reactions or digestions.
Loosen lids when thermally sterilising.

Bottles, wash, multi-lingual, safety venting, wide neck, LDPE



- Multi-lingual wide neck wash bottles featuring the patented DripLok® vapour venting valve that prevents solvent drips
- Indelibly printed in 4 languages (English, French, German and Spanish)
- Choose from five solvent labels all including the new Globally Harmonised Symbols, NFPA diamond, risk phrases and safety information
- Fitted with colour coded one-piece PP closures that maintain their shape in use
- Closures feature precision jet-tip aperture for fine control

Product Code	Capacity (ml)	Print label	Closure colour	Pack qty.
WGW531VTML	250	Acetone	Red	5
WGW535VTML	250	Ethanol	Orange	5
WGW534VTML	250	I.M.S	Yellow	5
WGW533VTML	250	Isopropanol	Blue	5
WGW532VTML	250	Methanol	Green	5
WGW537VTML	500	Acetone	Red	5
WGW541VTML	500	Ethanol	Orange	5
WGW540VTML	500	I.M.S	Yellow	5
WGW539VTML	500	Isopropanol	Blue	5
WGW538VTML	500	Methanol	Green	5

