



Aluminum Weighing/Drying Pan is ideal for general-purpose weighing, dispensing, storage, food processing, evaporating, or dust cover.

Additional applications include dispensing smaller quantities of liquids that may become usable after being exposed to the air for a period of time such as adhesives and epoxies. You can also use the dishes to hold small parts like hardware or small electronic components when manufacturing assemblies.

These dishes made from no oil residue and no vinyl coating to create contamination problems. The smooth bottom facilitates sample handling. Crimped walls increase rigidity and a large finger tab for sample numbers that can be etched prior to weighing for identification.

Capacity: 60ml, Dia.102 x Dep.8mm; 10boxes of50pcs per box

Capacity: 50ml; Size: Dia.120 x Dep.5mm. 10boexs of100pcs per box



20ml, Dia.45mm x Dep.15mm;

42ml, Dia.60mm x Dep.16 mm;

55ml, Dia.70 x Dep.16mm;

Wrapped in packs of 100; case of 1000 pieces.