




Nálevka Buchnerova  
 Funnel Buechner  
 Büchnertrichter  
 Embudo Buechner  
 Воронка Бюхнера



Označení Sign	D 1 mm	D 2 mm	D 3 mm	H 1 mm	H 2 mm	V/ml approx.	Filter paper diameter mm	Číselník Code JIPO
237/0	29	22	1	17	28	10	27	641 332 237 001
237/0a	42	33	1	20	35	25	40	641 332 237 011
237/1	48	38	1	24	43	35	45	641 332 237 101
237/2	62	45	1	30	64	70	55	641 332 237 201
237/3	77	58	1	35	64	120	70	641 332 237 301
237/4	97	70	2	40	71	240	90	641 332 237 401
237/5	116	95	2	49	83	400	110	641 332 237 501
237/6	130	110	2	52	85	600	125	641 332 237 601
237/7	156	130	2	59	96	1000	150	641 332 237 701
237/8	192	160	2	74	117	2000	185	641 332 237 801
237/9	248	220	2	90	117	4000	240	641 332 237 901
237/9a	296	260	2	105	130	7000	270	641 332 237 921
237/10	334	300	2	120	133	10000	320	641 332 237 991