



Combustion crucible - middle form



Sign	D (mm)	H (mm)	V (ml)	DIN
2/30	30	25	10	
2/35	35	28	12	DIN 12 904
2/40	40	32	20	DIN 12 904
2/45	45	36	30	DIN 12 904
2/50	50	40	45	DIN 12 904
2/60	60	48	80	DIN 12 904
2/70	70	56	120	DIN 12 904
2/80	80	62	200	

Combustion crucible - high form



Sign	D (mm)	H (mm)	V (ml)	DIN
3/30	30	38	15	DIN 12 904
3/35	35	44	26	DIN 12 904
3/40	40	50	35	DIN 12 904
3/45	45	56	50	DIN 12 904
3/50	50	62	72	DIN 12 904
3/60	60	75	130	DIN 12 904