

Sänger PRIMA Neotril, Neotril Gloves, Pack of 100

Description

Sänger PRIMA Neotril, Neotril Gloves, powder free, purple, Pack of 100

This unique polymer blend furnishes this glove with the comfortable fit of a latex glove and at the same time offers resistance of a nitrile/neoprene glove against many chemicals. The test list can be submitted upon request. Besides the medical area also specially suitable for work in the laboratory and in research. Additional treatment for easy donning. Lowest protein contents below 1.0 mg/g according to US-LEAP assys test. Slip resistant due to special surface treatment.

Features

- neoprene-nitrile-natural latex blend
- resistant to many chemicals
- powder free
- micro-textured surface
- very low protein contents
- easy donning, reduces sweating inside the gloves
- highly tear resistant
- superior tactile properties
- excellent fit
- rolled cuff edge
- length: 24 cm

Standards

- CE
- EN 455 part 1, 2 and 3
- EN 374, EN 420
- AQL 1.5
- EN ISO 9001 / 12.2000
- EN ISO 13485 / 2003
- FDA 1000-ml-water leak test
- tested for viral resistance and resistance against micro-organisms, suitability for use with food stuff

