

**W.A. Rubbermate Co., Ltd**  
**PRODUCT SPECIFICATION - NR-FT-53**

THE STANDARD OF LATEX AND COMPOUNDING MATERIAL SAFE FOR USE  
**5.3±0.2g 240MM POWDER FREE**  
**LATEX DISPOSABLE GLOVE - FINGER TEXTURED - NATURAL**

SKU Series:	NR-FT-53	Texture:	FINGER TEXTURED
Material:	LATEX	Texture (Cuff):	SMOOTH
Donning:	POWDER FREE	Texture (Palm):	SMOOTH
Length:	240MM	Texture (Finger):	TEXTURED
Weight Class:	5.3g±0.2g (Size M - Median)	Color:	NATURAL

Product Code	Size	Weight (g)	Single Wall Thickness (mm)		
			Cuff	Palm	Finger
NR-FT-53-S	SMALL	4.8 ± 0.2	0.10 ± 0.01	0.12 ± 0.01	0.14 ± 0.01
NR-FT-53-M	MEDIUM	5.3 ± 0.2	0.10 ± 0.01	0.12 ± 0.01	0.14 ± 0.01
NR-Ft-53-L	LARGE	5.8 ± 0.2	0.10 ± 0.01	0.12 ± 0.01	0.14 ± 0.01

<u>TEST STANDARDS</u>	LATEX	<u>NITRILE</u>
<u>Global</u>	ASTM D3578	ASTM D6319
<u>Europe</u>	EN455-1,2,3	EN455-1,2,3
<u>TENSILE STRENGTH</u>	LATEX	<u>NITRILE</u>
<u>Before Aging (mPa):</u>	Min. 18	Min. 14
<u>After Aging (mPa):</u>	Min. 14	Min. 14
<u>ELONGATION</u>	LATEX	<u>NITRILE</u>
<u>Before Aging (%):</u>	Min. 650	Min. 500
<u>After Aging (%):</u>	Min. 500	Min. 400
<u>POWDER RESIDUE</u>	LATEX	<u>NITRILE</u>
<u>Powdered (mg/dm<sup>2</sup>)</u>	Max. 10.0	Max. 10.0
<u>Powder Free (mg/glove)</u>	Max. 2.0	Max. 2.0
<u>PROTEIN RESIDUE</u>	LATEX	<u>NITRILE</u>
<u>Powdered (µg/dm<sup>2</sup>)</u>	Max. 200	0
<u>Powder Free (µg/dm<sup>2</sup>)</u>	Max. 50	0

<u>GRADE</u>	<u>PINHOLE AQL LEVEL</u>
Medical	AQL 1.5
Industrial	AQL 4.0
General Purpose	OFFLINE

<u>SIZES</u>	<u>PALM WIDTHS</u>
X-SMALL:	70 ± 5 mm
SMALL:	85 ± 5 mm
MEDIUM:	95 ± 5 mm
LARGE:	105 ± 5 mm
X-LARGE:	115 ± 5 mm
XX-LARGE:	125 ± 5 mm

Product Specifications provided are subject to change or availability. Please confirm with your account manager prior to finalizing your order.