



NEOTOUCH

A unique combination of allergy prevention and comfort.

KEY FEATURES AND BENEFITS

- The NeoTouch® gloves were the first neoprene single-use gloves for industrial usage.
- With a latex-free formulation, NeoTouch® gloves are appropriate for the prevention of Type I allergies. Because they're powder-free, they also limit the risk of dermatitis for the wearer.
- They provide excellent resistance to acids, bases and alcohols. Thanks to an exclusive formulation, NeoTouch® gloves are convincingly the most comfortable synthetic single-use gloves available today.
- NeoTouch® gloves have a polyurethane inner-coating that ensures easy donning. Their textured fingertips provide a secure grip in both dry and wet conditions. The beaded cuff ensures a secure hold, and their green colour makes them easy to identify.
- NeoTouch® gloves are inspected to an AQL of 1.5 for pinholes, which matches the legal requirements for medical gloves.

Product details

COLOUR

Bright green

EN SIZES

6.5 - 7,7.5 - 8,8.5 - 9,9.5 - 10

COATING MATERIAL

Neoprene

CUFF STYLE

Rolled beaded

GRIP DESIGN

Textured fingertips

THICKNESS (MM)

0.13

LENGTH (MM)

240

AQL

1.5

ackaging

100 gloves in a box; 10 boxes in a carton

up arrow

Standards and certifications

EN374:2003



EN374:2003



EN ISO 374-1:2016



KPT

EN ISO 374-5:2016



VIRUS

אלטרנטיבה לרגישים מגומי טבעי (לטקס)

מגינה מאלרגיות סוג 1

הגנה מעולה מחומצות, בסיסים ואלכוהוליים

מצויינת למעבדות

- כפפת נאופרן לעבודה עם כימיקלים וחומרים ציטוטוקסיים
- אורך 240 מ"מ, עובי 0.13 מ"מ
- עמידות גם מפני חומרים ציטוטוקסיים
- מתאים גם לאלרגיים ללטקס
- כולל תקן לכימיקלים, EN374 דרגת 1.5 AQL
- קיים במגוון מידות, מידה S עד מידה XL