semper**med®()** supreme

SUPREME QUALITY ALL ALONG THE LINE

EASY TO DON – patented synthetic inner coating SKIN-FRIENDLY – the Semperit post production process OPTIMUM SAFETY – ideal tactile properties



QUALITY ALL ALONG THE LINE

Sempermed[®] Supreme stands for quality and reliability; not without reason it is one of the most frequently used powder free surgical gloves in Europe.

EASY TO PUT ON

The lining does not form a smooth surface, but has a special net-like structure. This optimises the contact surface between the film and the skin, thus considerably reducing the friction resistance when putting on the glove. The structure leaves a pleasant, velvety feeling on the skin.

OPTIMISED COMFORT

Developed specifically for the requirements of the operating theatre, the fully anatomical shape with curved fingers corresponds with the natural, relaxed hand posture. The conical cuff is anatomically designed to fit the wrist, and the beading on the cuff provides a safe border to the gown.

TECHNICAL DATA

■ IDEAL TACTILE PROPERTIES

Due to its resistance, the wall thickness of the Sempermed[®] Supreme has been reduced to such a level that your tactile sensitivity is not in the least impaired, while at the same time providing optimal safety. The glove feels like a second skin.

SKIN-FRIENDLY

Sempermed has a specially developed, automatic finishing process in which the gloves are washed intensively a number of times. Skin-irritating substances resulting from the

production process and the natural latex protein concentration are thus reduced to the technically lowest level.

ТҮРЕ	powder free, with synthetic lining					
MATERIAL	natural rubber latex (NR)					
COLOUR	natural white					
GLOVE SHAPE	fully anatomical with rolled rim					
SIZE/OVERALL LENGTH as per EN 455-2	5½, 6 and 6½ / median 270mm 7, 7½ and 8 / median 280mm 8½ and 9 / median 285mm					
WALL THICKNESS measured in single layer in the palm area	0.19-0.23 mm					
BARRIER PERFORMANCE Tightness exceeding EN 455-1 (AQL 0.65 instead of AQL 1.5)	AQL 0.65					
FORCE AT BREAK as per EN 455-2	median 12N					
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years					
STERILISATION	radiation STERILE R					
PACKAGING	left and right glove turned inside out, in fibre-free inner pouch, ozone-tight, sealed in flat medical peel pack. in dispenser carton with sterilisation indicator: 50 pairs in transport carton with sterilisation indicator: 300 pairs					
LABELLING	EN 455-1/-2/-3/-4, EN 420, EN 421, EN ISO 374-1, EN 374-2, EN 556, EN 16523-1, EN 374-4, ISO EN 374-5, ISO 15223-1					
PURPOSE as set out in MD Directive 93/42/EEC as set out in PPE Regulation (EU) 2016/425	medical device class IIa single-use protective glove – PPE category III (protection against chemical substances for limited time)*					
ARTICLE NUMBERS	size 5½ size 6 size 6½ size 7	EU 822751521 822751601 822751621 822751701	INT 822751525 822751605 822751625 822751705	size 7½ size 8 size 8½ size 9	EU 822751721 822751801 822751821 822751901	INT 822751725 822751805 822751825 822751905

sterile surgical glove for single use.

CONTACT US!

Semperit Investments Asia Pte. Ltd. 8 Jurong Town Hall Road, #29–03/06 The JTC Summit, Singapore 609434, Singapore Tel.: +65 6274 4861 · Fax.: +65 6274 6977

Sempermed Europe GmbH Am Belvedere 10, A-1100 Vienna, Austria Tel.: +43 1 79 777-0 · Fax: +43 1 <u>79 777-630</u>

sempermed@semperitgroup.com · www.sempermed.com

*Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/chemical-resistance-brands **IMPORTANT NOTE:** The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper** storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used. In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** [©]Copyright 2022 Sempermed