



## GMC digital balance

1. Small in size and light in weight, fast weighing speed, high accuracy & stability
2. High quality and inexpensive, simple easy to operate, can replace plate frame balance
3. LCD display screen, more clear and bright
4. External CAL function
5. Various weighing units- gram、carat、 ounce etc
6. 6V 2.5A rechargeable battery
7. With RS232 interface
8. Perfect overload protection function
8. Optional with second display

Model	Capacity(g)	Readability(g)	Repeatability	Linearity	Operate Temp. (°C)	Pan Size(mm)	Housing Size (mm)	Response Time
GMC1002	1000	0.01	±0.01	±0.02	0-40	Φ180	230*270*70	≲2s
GMC2002	2000	0.01	±0.01	±0.02	0-40	Φ180	230*270*70	≲2s