



GMD digital balance

1. Small in size and light in weight, fast weighing speed, high accuracy & stability
2. High quality and inexpensive, simple easy to operate, can replace plate frame balance
3. LCD display screen, more clear and bright
4. External CAL function
5. Various weighing units- gram、carat、 ounce etc
6. 6V 2.5A rechargeable battery
7. With RS232 interface
8. Perfect overload protection function
9. Optional with second display

Product Detail

Model	Capacity (g)	Readability (g)	Repeatability	Linearity	Operate Temp. (°C)	Pan Size (mm)	Housing Size (mm)	Response Time
GMD1001	1000	0.1	±0.1	±0.2	0-40	Φ 180	230*270*70	≤2s
GMD2001	2000	0.1	±0.1	±0.2	0-40	Φ 180	230*270*70	≤2s
GMD3001	3000	0.1	±0.1	±0.2	0-40	Φ 180	230*270*70	≤2s
GMD5001	5000	0.1	±0.1	±0.2	0-40	Φ 180	230*270*70	≤2s
GMD6001	6000	0.1	±0.1	±0.2	0-40	Φ 180	230*270*70	≤2s