



**Cross**

A very stable general purpose stirrer



**Single Cross-Head**

Designed for use with tube-like containers yet very effective as general stirrers



**Double Cross-Head**

Similar to Single Cross-head but greater turbulence

Reference	D (mm) (D (in))	Reference	Ht (mm) x D (mm)	Reference	Ht (mm) x D (mm)
001.2401	10 ( 3/8 )	001.1110.1	8 x 10	001.1110	8 x 10
001.2402	20 ( 3/4 )	001.1114.1	12 x 14	001.1114	10 x 14
001.2403	25 ( 1 )	001.1117.1	13 x 17	001.1117	13 x 17
001.2404	30 ( 1 1/4 )	001.1125.1	15 x 25	001.1122	15 x 22
001.2405	38 ( 1 1/2 )	001.1140.1	17 x 40	001.1130	12 x 30
001.2406	50 ( 2 )	001.1160.1	17 x 60	001.1135	12 x 35
001.2407	60 ( 3 3/8 )			001.1140	14 x 40
				001.1160	15 x 60

NB: The dimensions of stirrer bars are nominal and approximate to the stated value by  $\pm 5\%$  for length and  $\pm 10\%$  for diameter