



cowie® PTFE Stirring Magnets

HOME

LABWARE
INDEX**Triangular**

Particularly effective for dissolving solids and mixing sediments because of the scraper like action on the base of the container

**Disc**

For use in tubes

Reference	L (mm) x D (mm)	Reference	Ht (mm) x D (mm)
001.412	12 x 6	001.709	6 x 9
001.420	20 x 8	001.710	6 x 10
001.425	25 x 8	001.720	10 x 20
001.425.14	25 x 14	001.730	12 x 30
001.435	35 x 9		
001.440	40 x 14		
001.450	50 x 12		
001.455	55 x 14		
001.480	80 x 18		

NB: The dimensions of stirrer bars are nominal and approximate to the stated value by $\pm 5\%$ for length and $\pm 10\%$ for diameter

cowie®